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### TESTING PROCEDURE (TP)

#### PH (LIQUID OR SLURRY SAMPLES)

ELECTRODES OR FOR ELECTRODE SOAKING BETWEEN USES.

1. 0.1M HCL: ADD 8.3 ML OF HCL (36%) TO WATER AND DILUTE TO 1 LITER WITH WATER. MIX AND LABEL CONTAINER "0.1 M HYDROCHLORIC ACID."
2. 0.1M NAOH: WEIGH 4.0 OF NAOH PELLETS INTO 1 LITER FLASK, ADD WATER TO DISSOLVE AND DILUTE TO VOLUME WITH WATER. MIX AND LABEL CONTAINER "0.1M SODIUM HYDROXIDE."
3. 0.1M KCL: WEIGH 7.45 G OF KCL INTO 1 LITER FLASK, ADD WATER TO DISSOLVE AND DILUTE TO VOLUME WITH WATER, MIX AND LABEL CONTAINER "0.1M POTASSIUM CHLORIDE."

#### SAMPLING

-----  
REPRESENTATIVE SUBSAMPLES FROM LOTS OR BATCHES SHOULD BE TAKEN. WHERE LARGE PARTICLES OR NON-HOMOGENEOUS FOOD SAMPLES ARE TAKEN. BLEND USING A CUISINART OR OTHER APPROPRIATE BLENDER/HOMOGENIZER.

#### PROCEDURE

##### A. STANDARDIZATION OF PH METER/ELECTRODE SYSTEM.

NOTE: THE PH METER/ELECTRODE SYSTEM MUST BE STANDARDIZED BEFORE EACH SERIES OF MEASUREMENTS AND RESTANDARDIZED AFTER 10 PH DETERMINATIONS (THIS INSURES THAT THE METER/ELECTRODE SYSTEM IS NOT DRIFTING). STANDARDIZATION STEPS GIVEN IN THIS TP ARE TYPICAL FOR MOST PH METERS THAT USE TWO-BUFFER STANDARDIZATION. WHEN THIS PROCEDURE DIFFERS FROM THE MANUFACTURER'S INSTRUCTION, THE LATTER SHOULD PREVAIL. IT IS ASSUMED THAT ALL MEASUREMENTS ARE MADE AT ROOM TEMPERATURE OR THAT BOTH SAMPLE ANALYSIS AND STANDARDIZATIONS ARE MADE AT THE SAME TEMPERATURE. VARIATIONS IN TEMPERATURES WILL CAUSE ERRORS IN PH MEASUREMENTS.

1. WARM UP PH METER FOR 1/2 HOUR (MOST METERS ARE LEFT IN STANDBY OR READY MODE WHEN NOT IN USE).
2. SELECT TWO STANDARD BUFFER SOLUTIONS AND POUR ABOUT 25 ML OF EACH INTO SEPARATE 50 ML BEAKERS. SELECT BUFFER SOLUTIONS WHICH WILL BRACKET THE ANTICIPATED PH OF THE SAMPLE. I.E., FOR SAMPLE PH'S LESS THAN 7.00, SELECT PH 4.00 AND PH 7.00 BUFFERS. FOR SAMPLE PH'S GREATER THAN 7.00, SELECT PH 7.00 AND PH 10.0 BUFFERS. USE ONLY FRESH BUFFER SOLUTIONS FROM BOTTLES.
3. SET TEMPERATURE COMPENSATOR CONTROL ON THE INSTRUMENT TO THE

TESTING PROCEDURE (TP)

## PH (LIQUID OR SLURRY SAMPLES)

OBSERVED TEMPERATURE OF THE BUFFER SOLUTION (CA 25 DEG.C.).

4. PLACE THE ELECTRODE(S) IN THE PH 7.00 BUFFER. GENTLY STIR THE SOLUTION BY SWIRLING AND THEN TURN METER SWITCH TO READ PH. ALLOW METER READING TO STABILIZE (UP TO 30 SECONDS).
5. SET THE METER TO PH 7.00 WITH THE "STANDARDIZE" OR "CALIBRATE" CONTROL.
6. RETURN METER SWITCH TO STANDBY POSITION. RINSE THE ELECTRODE WITH DISTILLED WATER AND BLOT DRY WITH A TISSUE. DO NOT RUB ELECTRODE WITH TISSUE. PLACE ELECTRODE IN THE PH 4.00 OR 10.00 BUFFER AS INDICATED IN STEP 2. GENTLY STIR SOLUTION AND TURN METER SWITCH TO READ PH. ALLOW METER READING TO STABILIZE (30 SECONDS).
7. ADJUST METER READING TO THE SECOND PH VALUE (PH 4.00 OR PH 10.00) WITH THE "SLOPE" ADJUSTMENT. NOTE: DO NOT READJUST THE CALIBRATE CONTROL FOR THIS SECOND PH BUFFER. IF STABLE READINGS CANNOT BE OBTAINED WITHIN 1 MINUTE, OR IF CHANGES OF MORE THAN 0.15 PH UNITS WERE REQUIRED AT EITHER THE PH 4.00 OR 10.0 ADJUSTMENTS, CONSIDER REJUVINATION OF THE GLASS ELECTRODE. SEE REMARK 3 FOR CORRECTIVE STEPS.
8. PLACE METER IN STANDBY OR READY MODE. DISCARD BUFFER SOLUTIONS WHEN FINISHED WITH STANDARDIZATION.

## PH DETERMINATION

1. TRANSFER 50 TO 100 ML OF SAMPLE TO A BEAKER FOR PH DETERMINATION.
2. INSERT STANDARDIZED ELECTRODE IN SAMPLE SOLUTION AND TURN METER SWITCH TO PH SETTING. READ AND RECORD THE PH FROM THE METER AFTER READING HAS STABILIZED.
3. RETURN METER SWITCH TO "STAND-BY" SETTING AND REMOVE ELECTRODE FOR SAMPLE. RINSE OFF ELECTRODE USING DISTILLED WATER. BLOT EXCESS WATER OFF WITH TISSUE BEFORE INSERTING ELECTRODE INTO NEXT SAMPLE.
4. WHEN ELECTRODES ARE NOT IN USE, STORE GLASS PH ELECTRODES IN PH 4.00 BUFFER. STORE SEPARATE REFERENCE ELECTRODES IN THEIR OWN FILLIGN SOLUTION. COMBINATION ELECTRODES SHOULD BE STORED IN A MIXTURE OF 50% PH 4.00 BUFFER AND 50% 0.1 M KCL

TESTING PROCEDURE (TP)

**PH (LIQUID OR SLURRY SAMPLES)**

SOLUTION. CHOOSE ELECTRODE STORAGE VESSEL THAT ALLOWS MINIMUM EVAPORATION.

CALCULATION

-----  
NONE

PRECISION

-----  
REPLICATE PH MEASUREMENTS SHOULD AGREE WITHIN +/- 0.005 PH UNITS.

REFERENCE

- 
1. J. ASSOC. OF ANAL. CHEM., 64, 2, 332 (1981)
  2. AOAC. OFFICIAL METHODS OF ANALYSES, 14TH EDITION, 1984, SECTION 32.010.
  3. DANGEROUS PROPERTIES OF INDUSTRIAL MATERIALS, 6TH EDITION, 1984, VAN NOSTRAND REINHOLD CO.

REMARKS

- 
1. ALL ELECTROLYTIC MEASUREMENTS REQUIRE A PAIR OF ELECTRODES - ONE TO SENSE THE ION OF INTEREST AND ONE FOR A REFERENCE OR KNOWN POTENTIAL. IN PH DETERMINATIONS, THE PH (OR GLASS) ELECTRODE IS PAIRED WITH ONE OF SEVERAL KINDS OF REFERENCE ELECTRODES. THE THREE MOST COMMON REFERENCE TYPES USED ARE THE CALOMEL, THE AG/AGCL AND THE ROSS ELECTRODE. OF THESE, THE ROSS REFERENCE ELECTRODE IS USED ONLY WITH A ROSS TYPE PH ELECTRODE. WHEN THE PH AND REFERENCE ELECTRODES ARE PLACED TOGETHER IN A SINGLE COMBINATION TYPE ELECTRODE BODY, ONLY THE AG/AGCL AND ROSS TYPE REFERENCE ELEMENTS ARE USED. THE COMBINATION TYPE ELECTRODE PROVIDES EASIER ACCESS INTO SMALL SAMPLE CONTAINERS AND SOME CONVENIENCE IN GENERAL HANDLING, CLEANING, ETC. BREAKAGE OR OTHER INACTIVATION OF ONLY ONE ELEMENT OF THE PAIR HOWEVER, ELIMINATES USE OF EITHER ELEMENT AND NECESSITATES REPLACEMENT OF THE ENTIRE COMBINATION ELECTRODE.
  2. ROSS IS A TRADEMARK NAME FOR ELECTRODES MANUFACTURED ONLY BY ORION RESEARCH INC. AND WHICH FEATURE A PATENTED DOUBLE-JUNCTION REFERENCE ELECTRODE AND PLATINUM CONDUCTING WIRES IN BOTH ELEMENTS OF THE ELECTRODE PAIR. THESE FEATURES PROVIDE SEVERAL ADVANTAGES OVER CONVENTIONAL CALOMEL OR AG/AGCL REFERENCE ELECTRODES. ROSS ELECTRODES FEATURE EXCELLENCE IN

TESTING PROCEDURE (TP)

PH (LIQUID OR SLURRY SAMPLES)

RESPONSE TIME, ELECTRODE STABILITY, RESISTANCE TO DRIFT, RESISTANCE TO TEMPERATURE EFFECTS ON PH MEASUREMENT AND RESISTANCE TO LONG TERM ELECTRODE FAILURE. WHILE MORE EXPENSIVE THAN CONVENTIONAL ELECTRODES, THE ADVANTAGES MENTIONED GENERALLY OUTWEIGH THEIR ADDED COST. SINCE NO SILVER OR MERCURIC CHLORIDE IS USED IN ROSS ELECTRODES, THE MOST COMMON CAUSE OF ELECTRODE FAILURE (PRECIPITATION OF AGCL OR Hg<sub>2</sub>CL<sub>2</sub> IN THE POROUS JUNCTION WITH RESULTING CLOGGING) HAS BEEN ELIMINATED. THE LONG-TERM LIFE OF ROSS ELECTRODES IS EXCELLENT.

ELECTRODE SYSTEMS THAT ARE SLUGGISH, ELECTRICALLY NOISY IN RESPONSE OR WHICH CONTINUALLY DRIFT USUALLY SUFFER FROM ONE OF THREE PROBLEMS:

- A. THE ELECTRODE IS CONTAMINATED FROM PRIOR SAMPLES ANALYZED. THIS CAN OCCUR WITH EITHER THE GLASS OR REFERENCE ELECTRODE (OR THE ENTIRE COMBINATION ELECTRODE).
- B. THE GLASS MEMBRANE ON THE PH ELECTRODE HAS CHANGED DUE TO SURFACE LEACHING. LEACHING OCCURS FROM BOTH NORMAL ELECTRODE AGING AND EXCESSIVE IMMERSION IN ALKALINE SOLUTIONS.
- C. THE REFERENCE ELECTRODE IS NOT PROVIDING A CONSTANT BACKGROUND POTENTIAL. WITH CONVENTIONAL ELECTRODES THIS IS USUALLY CAUSED BY A PARTIALLY PLUGGED POROUS FRIT OR JUNCTION.

CORRECTIONS OF THESE PROBLEMS CAN OFTEN BE MADE AS FOLLOWS:

- A. REMOVE GENERAL CONTAMINATION
  - 1) FOR PROTEIN OR INORGANIC DEPOSITS, WASH ELECTRODES IN 0.1 M HCL.
  - 2) FOR FAT OR GREASE DEPOSITS, WASH ELECTRODES WITH ACETONE OR ETHYL ETHER.

AFTER TREATMENT BY EITHER REAGENT, RINSE THE ELECTRODE THOROUGHLY IN DISTILLED WATER.

- B. RECONDITIONING A DAMAGED (LEACHED) GLASS PH ELECTRODE.

TESTING PROCEDURE (TP)

## PH (LIQUID OR SLURRY SAMPLES)

IMMERSE ELECTRODE TIP INTO 0.1 M HCL FOR ABOUT 1 MINUTE. RINSE WITH DISTILLED WATER THEN IMMERSE IN 0.1 M NAOH (0.1M KOH WILL WORK AS WELL) FOR ABOUT 1 MINUTE. CYCLE THE ELECTRODE THROUGH THESE SOLUTIONS SEVERAL TIMES ENDING WITH ELECTRODE IN ACID SOLUTION. AFTER THE FINAL DISTILLED WATER RINSE, SOAK ELECTRODE IN PH 7 BUFFER FOR 2 HOURS. RETEST FOR STABILITY AND ABILITY TO MEET SPAN REQUIREMENT (PROCEDURE STEP A-7). REPLACE ELECTRODE IF NECESSARY.

## C. UNBLOCK THE POROUS JUNCTION ON THE REFERENCE ELECTRODE.

## 1) CALOMEL REFERENCE ELECTRODES

DRAIN THE KCL SOLUTION FROM THE ELECTRODE CAVITY AND REPLACE WITH FRESH KCL SOLUTION. IMMERSE ELECTRODE IN WARM DISTILLED WATER AND PERIODICALLY APPLY AIR PRESSURE TO FILL HOLE TO RE-ESTABLISH ELECTROLYTE FLOW.

## 2) SILVER/SILVER CHLORIDE REFERENCE ELECTRODES

RINSE THE ELECTRODE CAVITY WITH DISTILLED WATER UNTIL ANY CRYSTALLIZATION OF ELECTROLYTE HAS BEEN DISSOLVED. DISCARD ALL RINSE WATER. REFILL CAVITY WITH FRESH 4 M KCL SATURATED WITH AGCL FILL SOLUTION AND APPLY AIR PRESSURE TO FILL HOLE TO ESTABLISH ELECTROLYTE FLOW.

IF ELECTRODE PERFORMANCE IS STILL UNSATISFACTORY, TRY THE FOLLOWING PROCEDURES IN THE ORDER LISTED AS REQUIRED:

- A) SOAK ELECTRODE OVERNIGHT IN 0.1 M KCL SOLUTION (NOTE: THIS SOLUTION IS MORE DILUTED THAN FILLING SOLUTION).
- B) BOIL THE JUNCTION IN DILUTE KCL SOLUTION FOR 5-10 MINUTES.
- C) CAREFULLY FILE OR SAND THE POROUS JUNCTION.
- D) DISCARD ELECTRODE

## CHEMICAL SAFETY INFORMATION

TESTING PROCEDURE (TP)

PH (LIQUID OR SLURRY SAMPLES)

- A. ETHYL ALCOHOL CAS 64-17-5 95% REF #3 PAGES 316-17

FLAMMABLE LIQUID, AVOID OPEN FLAMES, SPARKS OR HIGH HEAT SOURCES, USE IN HOOD TO AVOID VAPOR FUME BUILD UP.

NOTE: ETHANOL MAY BE DENATURED WITH METHANOL AND ISOPROPYL ALCOHOLS. DO NOT INGEST AS METHANOL AND ISOPROPYL ALCOHOL ARE TOXIC.

- B. HYDROCHLORIC ACID CAS 7657-01-0 REF #3 PAGE 1545

CORROSIVE LIQUID. AVOID SKIN CONTACT OR BREATHING VAPORS. USE IN A HOOD. USE PROTECTIVE EQUIPMENT TO PROTECT SKIN AND EYES SUCH AS GLOVES AND PROTECTIVE GOGGLES OR FACE SHIELDS AS WELL AS APRONS. WHEN DILUTING ADD ACID TO WATER NOT THE REVERSE. DO NOT MIX WITH ORGANIC COMPOUNDS. READ REFERENCE FOR POTENTIAL REACTIONS.

- C. SODIUM HYDROXIDE CAS 1310-73-2 REF #3 P. 2434, 2435

CORROSIVE MATERIAL. AVOID CONTACT WITH ANY BODY TISSUE. WEAR PROTECTIVE EQUIPMENT SUCH AS GLOVES, EYE PROTECTION AND APRONS TO PREVENT SKIN OR EYE CONTACT. WASH ANY EXPOSED AREAS WITH WATER. USE GOOD LABORATORY PRACTICES IN HANDLING THIS CHEMICAL.

- D. POTASSIUM CHLORIDE CAS 7447-40-7

SALT: NO PARTICULAR PROBLEM OR TOXICITY IN WORKING WITH THIS MATERIAL. DO NOT INGEST. USE GOOD LABORATORY PRACTICES.

ALTERNATE T.P.'S

CODE      DESCRIPTION

NONE

29.002/C/APPROVED  
October 11, 1990

LAST ISSUE CHANGE LOG  
TESTING PROCEDURE (TP)

PH (LIQUID OR SLURRY SAMPLES)

CHANGE LOG FROM LAST ISSUE:

CONVERSION TO THE NEW DOCUMENTATION SYSTEM.

CHANGE FROM LAST ISSUE

-----  
ITEMS CHANGED: TOTAL

REASON FOR CHANGE: TOTAL REVISION TO INCLUDE ELECTRODE CARE  
AND CLEANING, STANDARDIZATION AND  
CHEMICAL SAFETY INFORMATION.

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